NASAL WASH



• 28083 SINUS C SALINE NASAL SPRAY 20 ml - for children - box of 48 • 28084 SINUS A NASAL WASH - for adult - bottle 120 ml + 12 saline sachets Nasal wash helps to remove the sinus symptoms such as runny nose, nasal dryness and nasal congestion. By performing a nasal wash, you can remove particles such as pollen and dust in order to reduce the mucous membranes inflammation and make it

HUMI-RAINBOW

adjusted in direction (360° rotative head) and speed (adjustable knob)






• 28070 HUMI-RAINBOW HUMIDIFIER A cold-mist ultrasonic humidifier with

colours change • 28071 BALSAMIC

ESSENCE - pine 28072 BALSAMIC **ESSENCE** - citrus

28070



TECHNICAL SPECIFICATIONS

Power feeding: AC 110-240V / 50-60 Hz Max mist output: 280 ml/h Water Tank capacity: 2.8 l Suitable for room size: min 10, max 25 m² Max operating hours (non-stop): 10 hours Size (individual box): Ø 19xh 30 cm - Weight: 1.4 kg VAPINAL INHALER



• 28110 VAPINAL THERMAL WATER INHALER Warm jet mineral water inhaler ideal for home spring water inhalation treatment. It delivers hot humid steam through a nozzle working with the Venturi system and mixing the steam with the principles of thermal water. Equipped with stainless steel boiler working with resistance at 230V-50 Hz. Supplied with polycarbonate ampulla. Made in Italy.

• 28112 VAPINAL AMPULLA - spare part

TECHNICAL SPECIFICATIONS

Capacity: 800 ml Steam boiler in stainless steel

Security valve and thermostat Operating voltage: 230 V - 50 Hz - 600 Watt

Heating time: 5-8 min

Inhalation time (for 1 l): 12-15 min

ACCESSORIES FOR ANY BRAND OF NEBULIZER

ACCESSORIES FOR NEBULIZER

more resistant to allergies.

Gima offers different nebulizer bulbs with their line of accessories, suitable for any brand of nebulizer.



GIMA code	NEBULIZER ACCESSORIES	To be used with	
Standard accessories with Famiglia, Corsia, Eolo			
28139 28146	PVC connection tube - 1m Nebulizer bulb Facemask - adult* - replaces 28143	all nebulizers all nebulizers bulb 28139/28140	
28147 28157	Facemask - paediatry* - replaces 28144 Mouthpiece - nasal prong	bulb 28139/28140 bulb 28139/28140	
28138 28140	Nebulizer kit adult and pediatric facemask, pvc-tube, bulb, r Nebulizer bulb	all nebulizers prong, mouthpiece all nebulizers	

FASTERJET BULB with Valve System

A new designed Jet Nebulizer capable of delivering any medication in a very short time with a high amount of respirable particles. The Valve System makes therapy extremely efficient by reducing the amount of wasted medication while optimizing the output to the patient.

TECHNICAL SPECIFICATIONS FASTERJET

Capacity: 16 ml Neb rate: 0.60 ml/min Respirable volume: <5 µm 80%

Residual volume: 0.8 ml MMAD (Mass Median Aerodinamic Diameter): 2.7 µm

Fasterjet bulb: For delivering any medication in a very short

28163 + 28165



GIMA code	NEBULIZER ACCESSORIES	To be used with	
28163	Fasterjet bulb*	all nebulizers	
	with Valve System		
28164	Nebjet bulb **	all nebulizers	
	with Double Venturi Effect		
28165	Mouthpiece bi-valve	bulb 28163	
28166	Mouthpiece	bulb 28164,	
28167	Nasal Prong	bulb 28164, 28163	
28168	Facemask-adult***	bulb 28164, 28163	
28169	Facemask-pediatry***	bulb 28164, 28163	
28145	Filter - spare	28130/1/2/3/4/5	
* Standard with Happyneb II. Hospyneb ** Standard with Happyneb I. III. Betaneb			

*With elastic