

## GIMA SHARPS CONTAINERS

### POCKET



25780

### CS LINE



25783

25784



### PBS LINE



25793

25790



25795



25792

25794

25791

Double opening system

Temporary closure

Definitive closure



Wide opening

Indents for needles removal

Temporary closure: 1 click  
Definitive closure: 2 clicks

Indents for needles removal with safety fins

Temporary closure

Definitive closure



- pocket size
- ergonomic shape
- easy to use with just one hand
- double opening system: needles and insulin pen removal, entire syringe disposal
- maximum filling level indicator
- no emission of any harmful gas during the incineration process
- temporary and permanent closure
- possible side fastening

- simple and practical snap-type closure
- easy manoeuvre for temporary opening with just one hand
- wide opening for the disposal of voluminous material
- mechanism for the unhooking of insulin needles
- maximum filling level indicator
- no emission of any harmful gas during the incineration process

- ergonomic handle for the more spacious models
- anti-defluxion wings that prevent the outflow of the waste
- mechanism for unhooking insulin needles
- maximum filling level indicator
- no emission of any harmful gas during the incineration process

GIMA code	POCKET SHARPS CONTAINER	Minimum order	GIMA code	PBS SHARPS CONTAINERS	Minimum order
25780	Pocket 0.4 l - 12.3x6.2x16.5h cm	box of 70	25790	PBS 0.6 l - Ø 10x14.5h cm	box of 100
<b>CS LINE SHARPS CONTAINERS</b>			25791	PBS 1.5 l - Ø 14x16.7h cm	box of 100
25783	CS 3 l - 16.7x16.7x22.2h cm	box of 50	25792	PBS 2 l - Ø 14x21h cm	box of 75
25784	CS 4 l - 16.7x16.7x28.5h cm	box of 40	25793	PBS 5 l - Ø 23x19.5h cm	box of 55
<b>CERTIFICATIONS</b>			25794	PBS 7 l - Ø 24.8x19.5h cm	box of 40
FDA K071517			25795	PBS 12 l - Ø 24.8x31.34h cm	box of 25

## WASTE CONTAINERS - 30-50-60 l

### WASTE CONTAINERS - 3 sizes

- extremely solid and shock resistant at high and low temperature
- hermetic seal by peripheral adhesive on the lid (single or double)
- lateral handles
- excellent stacking possibilities

GIMA code	WASTE CONTAINERS			
	Tide	Lid	Max load kg	Size
25774	Container 30 l	single	15	40x30x39 cm h
25775	Container 30 l	double	15	40x30x39 cm h
25776	Container 50 l	single	20	40x30x56 cm h
25777	Container 50 l	double	20	40x30x56 cm h
25778	Container 60 l	single	25	40x30x65 cm h
25779	Container 60 l	double	25	40x30x65 cm h



Single lid with adhesive



Double lid with adhesive



## BD SHARPS CONTAINERS

### SHARPS CONTAINERS BECTON-DICKINSON

Space saving containers, autoclavable at 121°C, with safety cover. Conform to BS7320:1990, NF x 30-500, ADR UN 3291

- 25822 BD SHARPS CONTAINER 1.5 l - 16.5x16xh 20.5 cm
- 25823 BD SHARPS CONTAINER 3 l - 18.2x18.1xh 26 cm
- 25827 WALL SUPPORT for 25823
- 25828 BD SHARPS CONTAINER 7 l - 30.4x22.2xh 25.5 cm
- 25829 BD SHARPS CONTAINER 22.7 l - 32x22.7xh 44.5 cm



25829



25823



25828



25822

